

**AFFINITY HIS-TAG PURIFICATION
TECHNICAL SPECIFICATIONS
Empty Chromatography Spin Columns**

COLUMN MATERIAL	Polypropylene
FRIT MATERIAL	Polyethylene
FRIT PORE SIZE	35µm
CAPS	Two caps: Luer-lock & screw cap
CHEMICAL STABILITY	Stable in all commonly used reagents Autoclavable at 110°C
CATALOG NUMBER	P-300
CAPACITY	800 µL

For laboratory use only. Not for use in diagnostic or therapeutic procedures.

Gold Biotechnology

St. Louis, MO

Web: www.goldbio.com

Ph: (314) 890-8778

email: contactgoldbio86@goldbio.com