

## [rMuCSF2 Basic Information](#)

### **Amino Acid Sequence (124aa):**

APTRSPITVT RPWKHVEAIK EALNLLDDMP VTLNEEEVV SNEFSFKKLT CVQTRLKIFE QGLRGNFTKL  
KGALNMTASY YQTYCPPTPE TDCETQVTY ADFIDSLKTF LTDIPFECKK PVQK

### **Reconstitution:**

The sample should be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in a siliconized tube using PBS that contains a 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Reconstituted solutions are stable for up to one week at 2-8°C. Stock solutions should be aliquoted and stored at -80°C. Further dilutions should be made in appropriate buffered solutions containing BSA or serum.

### **Storage:**

This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage. The reconstituted sample can be apportioned into working aliquots and stored at -80 °C for up to 6 months. **Avoid repeated freeze/thaw cycles.**