

Stock Solution



TD-S Revision 2.0

Creation Date: 8/29/2014

Revision Date: 10/22/2019

2-10 mg/ml Cefsulodin Stock Solution

Instructions

1. Weigh 0.02-0.1 g of Cefsulodin (SCE-129) ([Cefsulodin, GoldBio Catalog # C-114](#), [CAS # 52152-93-9, mw. = 554.53 g/mol]).
2. Add 10 ml of sterile H₂O. Dissolve completely.
3. Prewet a 0.22 µm syringe filter by drawing through 5-10 ml of sterile H₂O and discard water.
4. Sterilize Cefsulodin Stock through the prepared 0.22 µm syringe filter.
5. Stock may be kept at -20°C for 1 year.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com