

Stock Solution



TD-S Revision 2.0

Creation Date: 10/15/2015

Revision Date: 4/12/2019

100mM Benzamidine HCl Stock Solution

Instructions

1. Weigh 174.63 mg of Benzamidine Hydrochloride Monohydrate ([Benzamidine Hydrochloride Monohydrate, GoldBio Catalog # B-050](#) [CAS 1670-14-0, mw. = 174.63 g/mol]) into a 15 ml tube.
2. Fill to a final volume of 10 mL with pure sterile H₂O.
3. Mix until completely dissolved.
4. Use immediately.

Note: Benzamidine is effective at a concentration of 0.5-5mM. Dilute the stock solution 1:100 to obtain a working concentration of 1mM.