

## dNTP Mix

Cat. No. D-900

GoldBio's 2'Deoxy nucleoside 5'-triphosphate mix is an aqueous solution containing dATP, dCTP, dGTP and dTTP, each at a final concentration of 10mM. This mixture has been tested for performance in a wide size-range of PCR templates with *Taq* and *Pfu* DNA polymerases.

The dNTP solution is ideal for PCR, cDNA synthesis, DNA sequencing and labeling procedures, and it has been performance tested with Q-PCR and long PCR (30 kb).

---

**Concentration:** 10mM

**Purity:** >99% by HPLC

**Quality:** Certified free of nucleases (DNase, RNase) and phosphatases.

**Storage:** -20 °C

---

Gold Biotechnology, Inc.  
1328 Ashby Rd.  
St. Louis, MO, 63132

**Phone:** 1-800-248-7609  
**Email:** [techsupport@goldbio.com](mailto:techsupport@goldbio.com)  
[www.goldbio.com](http://www.goldbio.com)