

Stock Solution



TD-S Revision 2.0

Creation Date: 4/8/2019

Revision Date: 4/8/2019

5mM (200x) DMNPE-Caged Luciferin Stock Solution

Instructions

1. Dissolve 2.45 mg of DMNPE-Caged Luciferin ([GoldBio Catalog # L-130](#)) in 1 ml of DMSO for a 5mM (200x) stock solution.
2. DMNPE-Caged Luciferin can be added to cell cultures at a final concentration of 25 μ M or at a concentration specified by protocol.
3. Store stock solution at -20°C and protect from light.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com